

APPENDIX

Contents of Chapter.

- [The Pascal Headers](#)
- ["Standardless" Analysis](#)
- [User Interfaces](#)
- [Make File of n Noise Spectra](#)
- [Gaussian and Lorentzian Distributions](#)
- [Measure Beam Energy](#)
- [Hardware Interfacing](#)
- [F Keys and Keyboard Shortcuts](#)
- [Detector Window Transmission](#)
- [Background Filtering](#)

[Table of Contents](#)

APPENDIX I. The Pascal Headers

CONST

Maximum_Channels = 8192;

MaxNumPks = 50; { Maximum peaks in a FIT, or window }

TYPE

Spectrum_counts = ARRAY [1..Maximum_Channels] OF Real;

Experiment Header;

{ The following record, Expt_InfoRec, comprises the EXPERIMENT header. }

Expt_InfoRec = RECORD

LastSpect : Integer;

FirstSpec : Integer; { first Spectrum# }

Specimen_ID : String[50];;

MCA_Filename : String[25];;

Specimen_Comment_Field : Str255;

Was_PassWord : String[25];

RefFile: Boolean;

Analyst : String[50];

Detector : [Detector Rec](#);

Azimuth : real; { degrees }

Elevation : real; { degrees }

Detector_Area : real; { cm sq }

Detector_Thickness : real; { mm }

C_Thickness : real; { um }

Diamond : real; { spare }

Mylar : real; { um }

BN_Thickness : real; { um }

SiN_Thickness : real; { um }

Ice_Thickness : real; { um, as pure oxygen }

Au_Thickness : real; { um }

Al_Thickness : real; { um }

Be_Thickness : real; { um }

Si_Thickness : real; { um }

Moxtek : real; { um }

Paralene: real; { um }

WDS_Resolution: real; { eV }

dE : real;

Si_Resolution : real; { eV at Mn K alpha }

Energy_Intercept : real; { y intercept in eV }

```
Energy_Slope : real; { }  
Number_of_Channels : integer;  
kV : real;  
Detector_Tilt : real; { Beta }  
Quantum : real;  
Spare2 : Boolean;  
Spare3 : integer;  
END;
```

```
Detector_Rec = RECORD  
Spec : (EDS,WDS,EELS);  
ID : integer;  
END;
```

{The SPECTRUM header is made up by combining the following records.}

Spectrum Header;

```
Element_InfoRec = RECORD  
Atomic_number : Integer;  
spare1 : real;  
Weight_Fraction : real;  
spare2 : real;  
Valence : real;  
END;
```

```
Plot_InfoRec = RECORD
```

```
Plot_Connected : Integer;  
Plot_Symbol : Integer;  
Spectrum_Color : RGBColor; {Red, Green, Blue }  
END;
```

```
Acq_InfoRec = RECORD
```

```
Probed_Area : real;  
X_Position : real; { um }  
Y_Position : real; { um }  
Spare1 : real;  
FirstChannel : integer;  
LastChannel : integer;  
Begin_Faraday : real; { nA }  
End_Faraday : real; { nA }  
Begin_Time : Longint;  
FirstValue : real;  
EndValue : real;  
spare2 : integer;  
Real_Time : real; { The time on the wall }  
Live_Time : real;  
Slow_Channel_Counts : Longint; { Represents total out counts }  
Medium_Channel_Counts : Longint; { If UTW. Approx. input counts below 1 keV }  
Fast_Channel_Counts : Longint; { Represents total input counts above 1 keV }  
RequestedLiveTime : LongInt;  
ActualLiveTime : LongInt;  
Acquiring : Boolean;  
LLD : Integer; { Acquisition setup dialog box }  
Offset : Integer; { Acquisition setup dialog box }  
PulseProcessorType : Integer;  
PulseProcessorSetting : Integer;  
END;
```

Spectrum_InfoRec = RECORD

```

Spectrum_Type : String[4];
Spectrum_Comment_Field : Str255;
Spectrum_Number : Integer;
Spectrum_Class : String[25];
Theoretically_Generated : Boolean;
This_is_a_Standard : Boolean;
BkgSubtracted : Boolean;
Maximum_Counts : real;
Minimum_Counts : real;
X_Tilt : real; { degrees }
Y_Tilt : real; { degrees }
Take_Off_Angle : real; { degrees }
Spec_Detector_Distance : real; { mm }
Spare : real;
Specimen_Thickness : real; { cm }
Specimen_Density : real; { g/cm sq. }
Number_of_Elements : Integer;
Element_Info : array[1..15] of Element_InfoRec;
extra_space : array[1..157] of real;
WDS_in_eV : Boolean;
bool2 : Boolean;
Average_Z : real;
Spare1 : real;
Spare2 : Boolean;
Spare3 : Integer;
END;
```

Element_InfoRec = RECORD

```

Atomic_number : Integer;
Spare1 : Integer;
Weight_Fraction : real;
Spare2 : real;
Valence : real;
END;
```

{In the next two RECORD definitions, we combine the above records into the final Work_Spectrum RECORD. The "Results" spectrum, as well as the "scratch" spectra 1-8, are copies of the Work_Spectrum RECORD.}

Spectrum_Structure = RECORD

```

Spectrum_Info : Spectrum\_InfoRec;
Acq_Info : Acq\_InfoRec;
END;
```

Work_Spectrum = RECORD

```

Expt_Info : Expt\_InfoRec;
Plot_Info : Plot\_InfoRec;
SpectrumStuff : Spectrum_Structure;
S : Spectrum_counts; { ARRAY [1..Maximum_Channels] OF Real }
END;
```

A DTSA file of spectra will contain one Expt_InfoRec RECORD followed by as many Work_Spectrum RECORDs as there are "spectra".

[Top](#)

[Table of Contents](#)